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住所：四川省成都市天府新区新兴街道天工大道916号A2栋3~4层
Residence: 3rd-4th floor, Building A2, Tiangong Road No. 916, Xinxing street, Chengdu, Sichuan Province People's Republic of China

生产地址：四川省成都市天府新区新兴街道天工大道916号A2栋3~4层
Manufacturing address: 3rd-4th floor, Building A2, Tiangong Road No. 916, Xinxing street, Chengdu, Sichuan Province People's Republic of China

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C013-1 1mLx2

Explanation of Symbols

	Sterilized using aseptic processing techniques
	Batch code
	Date of manufacture
	Use-by date
	Temperature limitation
	Keep away from sunlight

Precautions and Warmings

- Before using the product for embryo culture, the user should read and understand the operation and precautions, and master the correct procedures.
- This product cannot be used for injection.
- This product is sterilized by filtration and is not suitable for re-sterilization
- Do not use this product if:
 - The product package is damaged or the seal has been destroyed.
 - The product exceeds the shelf life.
 - The product has been discolored, or become turbid, or granular.
- This product contains gentamicin sulfate. Appropriate preventive measures should be taken to ensure that patients are not allergic to the antibiotic.
- In order to avoid contamination, please use aseptic technology during operation.
- After taken out, the product cannot be poured back to the vial for further storage and use, and

the unused product after balancing should be discarded.

8. If the semen quality is good, do not use the product to avoid hyperactive movement.

9. This product has no effect on dead sperm.

10. The components of the product have no loss in the specified storage condition and shelf life.

11. This product contains human serum albumin derived from human blood. The human serum albumin in this product has been tested to be negative for HIV-Ab, TP, and Anti-HCV, and HBsAg, HBsAb, HBeAg, HBeAb, and HBcAb. All blood products should be considered potentially infectious. Appropriate precautions should be taken to ensure that the patients do not have an allergic reaction to the human serum albumin.

Intended Use

The product is used for in vitro treatment of sperm specimens obtained from testis and epididymis of those with severe oligozoospermia and asthenospermia, to improve sperm motility.

Components

L-Proline	L-Leucine
L-Histidine	L-Isoleucine
L-Phenylalanine	L-Threonine
L-Asparagine	L-Lysine monohydrochloride
L-Aspartic acid	L-Valine
L-Glutamic acid	L-Serine
L-Arginine	L-Methionine
L-Tryptophan	L-Cystine
L-Tyrosine	Sodium chloride

Potassium chloride	Calcium chloride
Glycine	Magnesium sulfate
Sodium bicarbonate	Potassium phosphate monobasic
Glucose	N2-L-Alanyl-L-glutamine
Sodium pyruvate	Human Serum Albumin
Sodium lactate	Gentamicin Sulfate
Pentoxifylline	

Storage and Transportation

Store dark at 2°C to 8°C, is stable until the expiry date shown on the labels.

Product can be used within 7 days from first opening when stored dark at 2°C to 8°C.

Store dark and closed at 2°C to 8°C during transportation.

Discard unused medium within 7 days of opening.

Performance

- pH(37°C): 7.2~7.6.
- Osmolality (mOsmol/kg) : 260~310.
- Endotoxin (LAL): <0.1 EU/mL.
- Sterility Test (bacterial and fungal screen, SAL 10⁻³): Pass.
- Human sperm survival assay: Pass.

Contraindications

Those who are allergic to human serum albumin and gentamicin sulfate are prohibited.

Operation

Intrauterine insemination (IUI)

According to the standard procedure, after

preparing the semen for artificial insemination, add SpermAct according to one tenth of the semen volume for 10-30 min, then conduct IUI.

Intracytoplasmic sperm injection(ICSI) procedure is for patients with MESA and oligoasthenospermia. In the final culture medium, add SpermAct one tenth of culture medium, culture in 37°C constant temperature environment for 10-30 min, find the moving sperm, and then put it into the sperm micro working fluid according to the normal ICSI procedure.

Testicular sperm extraction(TESE)

Culture medium is prepared according to the standard TESE procedure. In the final culture medium, add SpermAct one tenth of culture medium, culture in 37°C constant temperature environment for 10-30 min, find the moving sperm, and then put it into the sperm micro working fluid according to the normal ICSI procedure.

Specimen Freezing and Thawing

The frozen sperm is prepared according to the standard process. According to one tenth of the volume ratio, add SpermAct to the sperm washing liquid.

After washing twice, it can be used.